


Hollow Fiber Filter Modules Standard Product Number Key

K **06** - **E** **100** - **05** - **N**

① ② ③ ④ ⑤ ⑥

Module Family ①

C	MicroKros MLL x FL
D	MidiKros FL x FL
T	MidiKros TC 1/2" TC x FL
S	MiniKros Sampler 3/4" TC x 3/4" TC
N	MiniKros 1.5" TC x 3/4" TC
K	KrosFlo® 3" TC x 1.5" TC
X	KrosFlo® MAX 6" TC x 1.5" TC




Membrane Type ③

E	mPES (Modified Polyethersulfone)
S	PS (Polysulfone)
P	PES (Polyethersulfone)
M	ME (Mixed Cellulose Ester)

MWCO Rating ④

001	1 kD	UF
003	3 kD	
005	5 kD	
010	10 kD	
030	30 kD	
050	50 kD	MF
070	70 kD	
100	100 kD	
300	300 kD	
500	500 kD	
750	750 kD	
05U	0.05 µm	
10U	0.10 µm	
20U	0.20 µm	
45U	0.45 µm	
50U	0.50 µm	
65U	0.65 µm	



Packaging ⑥

N	Normal
S	Sterile
P	Pre-Wet