



Pneumatic Scale
Angelus

CARR Centritech
Separation Systems

Powerfuge® P6



Features

- Drier solids and clear liquids
- Sealed automatic operation
- Full CIP/SIP designs
- cGMP construction
- Designs available for ASME-rated boundary, hazardous locations and aseptic operations
- State-of-the-art manufacturing facility located in Clearwater, Florida USA

Powerfuge® P6

PneumaticScaleAngelus' unique Powerfuge® P6 System provides continuous-flow two-phase separation of solids and liquids at rates up to 60 liters per hour. Using an innovative design, the feed is introduced through a stationary pipe and accelerated to full rotational speed before entering the titanium-alloy bowl. In this system, centrifugal forces as great as 20,000g permit separation of even sub-micron particles. Compressed solids are periodically removed with a fully-automated scraping cycle. After the solids have been discharged, automated CIP/SIP cycles can be performed via the PLC-controlled operating sequence. The P6 is designed for ease of operation, with sealed feed and drain lines allowing for aerosol-free operation. Similarly, no special tools are required to disassemble the system. Integrated turnkey systems incorporating feed delivery options, CIP and SIP skids, temperature control, and closed-loop centrate pumping under vacuum are available. The Powerfuge® P6 system meets cGMP and CE requirements and has been validated to meet chemical and pharmaceutical industry protocols. The P6 is engineered to optimize ease of operation while maximizing throughput; each step in the process can be readily automated so that time spent on routine maintenance can be minimized.

For more information, please contact us at:

CellD

20bis, rue du Chapitre

F-30150 ROQUEMAURE

Tel : +33 (0)4 66 82 82 60

Fax : +33 (0)4 66 90 21 10

Email : contact-info@celld.com www.celld.com

