



Centrifuge Product Line Overview

PneumaticScaleAngelus features four scalable models of automated centrifuges to meet separation needs



Powerfuge®

- High G
- Pelleting
- Classification & Clarification
- Submicron particles
- 1,000-20,000 xg, 10-1,700 LPH



Viafuge®

- Low Shear
- Viable, mammalian, insect & plant cell harvesting
- Clarification
- 500-10,000 xg, 30-2,400 LPH



Centritech®

- Single-Use
- Sterile
- Disposable
- Perfusion
- Viable, mammalian, insect & plant cell perfusion



UniFuge®

- Single-Use
- Sterile
- Disposable
- Viable, mammalian, insect & plant cell harvesting
- 300-4,000 xg, 30-1,800 LPH



CellD

20bis, rue du Chapitre

F-30150 ROQUEMAURE

Tel : +33 (0)4 66 82 82 60

Fax : +33 (0)4 66 90 21 10

Email : contact-info@celld.com www.celld.com

For more information, please contact us at:

